



9ED'BUa Y.

Œ ǻ} dŮ^&[{ { ^} â^â

GYfj JW'HmdY.

UÛ

7 ca dUbmBUa Y.

.....

7 i glca Yf'BUa Y#ŋhY.

.....

9ED': JYbUa Y.

ÒÙÈĞİ FÈ~]

9ED'Di V'jg\ '8 UHY.

R' } ^/GEF

Df]bh8 UHY.

Œ *~• ǻ ĞĖGHĖK HKİ ÁŮ

HUV`Y`cZ7cbhYbŋ

GYWŋcb	DU`Y
Ù&[]^Áŋ áÁŰ`!][•^	H
Ô! [••ŠæÄÖ[{] æŋ &^	FF
œÖÒÄ^ ä`!`^ÁŰ] ä }•	FÍ
œÖÒÄ^ ä`!`^ÁŰ^ÄÖæ^•	Fİ
ÒÙUÛ	FJ
Û^][!áŋ áÄ^ ä`!`^ÁŰ] ä }•	Gì
Ù^ &c^áÁŰŋ }æ`!`^ÁŰ] ä }•	Gì
Ô`•ä { ^!Áŋ]![:çæŋ	GJ
Š^* æÄ[äŰ^	H€

[illegible]

I bXYghUbX]b[`k Y`HYghGdYWZVU]cb`GYW]cb`b`HUVi`U`F`Yj`Yk`8cW`a`Yb]g

[illegible][illegible]

Q Áááá | æ ÁÚÁÁ | | çæÁ | & { ^ } ð @ Á • ð ^ d [á • Á æ æ ^ c ! • ð æ á Á æ Á æ Á | æ ^ á Á Á @ Á || , á *
| ^ | ^ • ^ æ æ Á | { æ È

Di a d': `ck '5 WW fUWn

 \emptyset, \mathbb{A}

OŒ æ̂ cāæÁœ ā̂} óæg åÁ [] Æœ ā̂} óĀ ~ { }

K	Únd [ă d	€€ { ŠD ă	ŠK €€	PK í €€
	Únd [ă d	€€ { ŠD ă		
	Šă ăK	ŁM í €€ Ā	ŠK í €€	PK F€€

œ ă ʰ} dŭ ʰ&[{ { ʰ} ă ʰă
ò • è € Ğ Ĥ f è ˘]

$$U_{\alpha}^* \wedge A \wedge \neg A \in$$

OE *v • d Ā Ē Ĝ Ğ H Á K H K Î Á Ú

ÔÖÜK

Ô`•q{ ^!ÁææÁ`•c{

Þ^ç [\ \ ÁÖÖK

ÖÖÖÁ/Á•ç\^âÁ} Á/Áç [\ \ Á[â^Á ã@ Á@ Áæ[\æ \ ^ ÁÖÞ/Á ÷æ d` &c \^

Š[&æ/ÖÖK

ÖÖÖÁ•æ•Á} Á/Á/á^] ^} á^} c^c^} æÁ/á^Á} } ^&c^áÁ/ Á@ ÁÖÞ/Á ÷æ d` &c \^

Öc^} ææ^Á^ ^c@ áK

ÖÖÖÁ/Á•ç\^âÁ} Á/[{ ÁÁÜÖÁ[cÁ} } ^&c^áÁ/ Á@ Áæ[\æ \ ^ Á^ç [\ \

Üæ åæá/ÖÜÜK

ÖÜÜÁ/ ^É^÷ ^áÁ/ ÁÉ ð^} c^æ/ &` á^/Á^æ |c^••É^d [ð •É} áÁ ã

Ö[] -ä` \^á/ÖÜÜK

ÖÜÜÁ/æ/æ Á^} &æ *^á/æ^âÁ} Á•^/Á^` ä{ ^} •

Xææ & ÉÜæ *^K

V@Áæ *^Á Á} ^\ææ * Áæ ^•Á/ Á/Á•c^ ææ ^c^Á ã@ ÁÖÖÁç[, É{] ^\æ \^É} áÁ[Á} Dæ ^\^Á^÷ ^áÁ/ á/Á•c^á @} Á@ c^c^ æ Á/Á^ç^[[\^áÁ/ áÁ^c^]

Xææ & ÉY K

Öæ *^Á/ Á/Á•c^ ææ ^c^Á/æ/Á ã@ Á@ ÁÉ ð^} c^ææ & Áæ *^

Xææ & ÉÜK

Öæ *^Á/ Á/Á•c^ ææ ^c^Á/æ/Á` •æ^Á@ ÁÉ ð^} c^ææ & Áæ *^

Xææ & ÉÖY K

Ö`æÁ^[[c^ *Á ã Á ã@ ÁÉ ð^} c^`]] [\c^áÁæ *^

Xææ & ÉÖUK

Ö`æÁ^[[c^ *Á ã Á` •æ^Á@ ÁÉ ð^} c^`]] [\c^áÁæ *^

U] ç} æÁ^•c

ç Á[] ç} æÁæáæ} æ/Á•c^æ/Á [c &` á^á/ Á@ Áæ åæá/ÖÜÜ

Öæáæ} æÁ^d [ð c

Öæáæ} æÁ ^æ` \^ { ^} c^æ ^Á/ Á/Á•c^ ææ ^c^Áç[, É{] ^\æ \^É} áÁ[Á} D

Ei U]ZWU]cb`@Z'7mWYGHU Yg

WÜÜK

W^!ÄÜ^~ ã^ { ^} oÄ] ^&ãBæÄ } Á@æÄ [&~ { ^} oÄ ç) â^âÄ •d~ { ^} oÄ •^

ÖÜK

Ö^•ã } Ä~ æãBæÄ } Á@æÄ [&~ { ^} oÄ ~ ææãæ Ä -Ä •d~ { ^} oÄ !Ä ç) â^âÄ •^

ŬK

Ŭ•æ||æÄ } Ä~ æãBæÄ } Á@æÄ [&~ { ^} oÄ @Ä •æ||æÄ } Ä -Ä @Ä •d~ { ^} c

UÛK

U] ^!æÄ } æÄ~ æãBæÄ } Á@æÄ •oÄ @Ä æ *^Ä -Ä •^

ÚÛK

Ú^!- { { æ &Ä~ æãBæÄ } Á@æÄ •oÄ @Ä •ç { Ä } â^!Ä { } âæÄ } •Ä -Ä •^

ÜÛK

Ü^] æÄ~ æãBæÄ } Á@æÄ •oÄ @Ä •d~ { ^} oÄ ç!Ä ç Ä •d~ { ^} oÄ] æÄ

ÚT K

Ú|æ } ^âÄ æ ç) æ &Ä !Ä &@ã |âÄ !^ç^} çÄ^Ä æ ç) æ &^

FYZfYbW'Gci fWVg

œJVT K

œ ^!ææ Ä [&ç Ä !Ä^•ç *Ä æ âÄ æ^!æç

øøœK

Wj æ^âÄæ^Ä [[âÄ æ âÄi~ *Ä { æ æ çæÄ } }

ōœ Ūt K

ō [[âÄ ç { æ^âÄ æ ~ æç !æ *Ä !æææ^

ŬUK

Ŭç!} æÄ } æÄ!~ æ ä æÄ } Ä !Ä ç ææâä æÄ }

ŬUJ€€FK

ŬUÄ ç ææâÄ !Ä~ ææ Ä æ æ^ { ^} oÄ •ç { •

ÚŬŬK

Úœ { æ^ æçŬ •^ &ç } Ō [[^!æÄ } Ä &@ { ^

WÙÚK


Wj ã^âÁÛæ^•Áú@æ{ æð[]^ãæ



Agilent CrossLab Compliance Services

Agilent CrossLab Compliance Services are designed to seamlessly integrate with traditional quality systems used by firms and recognized by regulatory agencies worldwide. Analytical instruments must be suitable for their intended use. This requirement is good science in all laboratories and a regulatory requirement in pharma and biopharma laboratories. A life-cycle process for documenting and testing the suitability of laboratory instruments should be followed and Agilent recommends the life cycle framework defined in USP General Chapter <1058> on Analytical Instrument Qualification (AIQ). USP <1058> defines the governing framework and requirements that need to be satisfied, but the laboratory is responsible for how they satisfy these requirements.

- The United States Pharmacopeia (USP) is the only major pharmacopeia with a general chapter dedicated to analytical instrument qualification, making <1058> an important global regulatory reference. The information is provided in a scientific, risk-based approach to analytical instrument qualification (AIQ). However, the life-cycle framework contained within USP <1058> is not prescriptive in its implementation, making the embedded scientific and risk-based principles flexible and universally applicable.
- The scientific process followed by CrossLab uses the Agilent's Automated Compliance Engine (ACE) to deliver paperless electronic qualification. The life-cycle stages Agilent perform are highlighted in the life-cycle diagram below. As part of this life-cycle, Agilent can configure the qualification tests performed to align with user requirements.



USP <1058> AIQ Framework

NOTE: RQ services, described later in this document, can be added to standard qualification services.

ACE Workflow and Equipment Qualification Plans (EQPs)

Overview

Within the ACE workflow, the qualification tests, setpoints, and limits are defined in an EQP that can be configured to ensure that testing satisfies user requirements. When the qualification work is complete, an Equipment Qualification Report (EQR) is issued. The electronic workflow used within ACE has significant data integrity advantages over traditional paper or Excel-based qualification protocols, as validated calculations can be performed directly using electronic data such as chromatograms and metrology test values. Several of the instrument life-cycle stages are the responsibility of the laboratory, Agilent can provide compliance consultancy services and documentation which can help customers satisfy these requirements. These additional services are not included in our typical qualification offering.



High-level ACE Qualification Workflow

Standard and User-defined Limits

(Hardware qualifications only)

EQPs are available for download and approval as standard documents with Agilent recommended tests, setpoints, and limits, or they can be electronically configured by approved personnel to align with user requirements and intended range of use requirements. The degree of configuration depends on the analytical technology, but most EQPs can be configured to some degree, and one feature that can typically be changed is test limits.

EQPs are designed to be configurable (dependent on the analytical technology and standard requirements), but including additional tests or setpoints can impact the qualification time and associated cost. If a test limit is changed, ACE includes the capability to report results against the Agilent approved limit and any customer required limits (that is, both can be reported simultaneously).

If a user-defined test limit is more stringent than an Agilent recommended limit, Agilent makes no guarantee or obligation regarding the instrument passing the tighter test specification requirements. It is important to appreciate that tests performed under conditions of use (that is, to satisfy pharmaceutical monograph and application requirements) can have different limits than those defined in the QQ. It is the continuum of the combined QQ, PQ, and any point of use testing performed each time the instrument is used that together satisfy regulatory requirements.

User Requirements Specification (URS)

The purpose of user requirements is to document the intended use of the instrument within the life-cycle process and quality management system (QMS) being followed. Therefore, the URS is a customer / laboratory responsibility. Defining user requirements is often used to guide the customer in instrument selection and is stated as the first activity that should be followed in <1058>. The URS is important for two main reasons.

- It is a regulatory requirement for FDA and EU GMP that the intended use of the instrument and any software must be specified.
- Investment protection perspective means getting the right instrument for the right job.

Qualification protocols should test the instrument against any limits or specifications listed in the URS, which should document the intended range of use. Depending on the instrument complexity and how it is classified, a separate URS document may not be needed, but the URS requirements of the <1058> framework must be satisfied. A separate URS is almost always recommended for computerized systems.

An instrument performance specification is a product of the instrument development process by the supplier. It typically documents the performance the instrument can achieve. The URS should be based on intended use of the instrument and not the instrument specification. Additionally, if the intended use of a system changes, this may trigger a need to review the URS and associated qualification testing (for example, to ensure range of use is tested if used with a new analytical procedure).

Agilent offers compliance consultation services and documentation that can help customers address URS requirements.

Design Qualification (DQ)

The main function of the DQ stage of the laboratory instrument life-cycle process is to document why the selected instrument is suitable. Typically, this includes consideration of the instrument specification, how the instrument will be qualified, and the QMS followed by the instrument manufacturer. All together, these confirm that instrument performance is capable of satisfying user requirements. Depending on laboratory instrument life-cycle policy or SOPs being followed, instrument requirements and the relationship between the URS and DQ stages may vary – but as long as the <1058> framework principles are satisfied, this is not a problem, as it is left to each laboratory to justify and document its specific approaches.

The responsibility for satisfying DQ requirements primarily lies with the laboratory, with support from the supplier.

Agilent's approach to satisfying DQ requirements of USP <1058> includes the following.

- All Agilent hardware and software laboratory products, including the ACE software used to deliver qualification services, are designed, manufactured, and tested according to Agilent internal quality life-cycle development procedures.
- Certificates of Agilent testing, validation, and conformance to standards are provided with new Agilent instruments and similar certification can be provided for ACE software.
- Agilent is capable of installation, support, preventive maintenance, on-going qualification, and re-qualification after repair and user training worldwide.

Agilent offers a compliance consultation service that can help customers with DQ documentation.

Agilent offers service contracts to repair and requalify an instrument during the period between scheduled annual OQs.

The level of retesting is prescribed in the RQ section of ACE: a form is displayed for the operator showing all types of repairs possible and the retesting required. Part of an example form for an LC system is shown below.

Re-Qualification After Repair		
Pump Strategies		
Repair/Replace Strategy	Modules	OQ Testing
Internal pump head parts, active inlet valve (or AIV cartridge), (parts of) check valves, reference valves, inlet manifold or pump drive, or taking pump head apart to clean (versus repair)	Any pump	Flow Accuracy & Precision
Pulse damper, pressure transducer	Any pump	Flow Accuracy & Precision
Multi-channel gradient valve	Quaternary	Flow Accuracy & Precision Gradient Composition

The full list of RQ repair and retest guidance is available for customer review.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2023
Published in USA



Overview

Agilent recommends use of **Network ACE** for CrossLab qualification services that are enabled using the Agilent Automated Compliance Engine (ACE) software. Network ACE and Local ACE both access data directly (default methods) and are considered equivalent from a data integrity and data traceability perspective (see below). To provide additional flexibility in qualification service delivery, an alternative method is also available that accesses data indirectly. Use of the alternative method requires customer pre-approval using this form.

Available Methods

Method	Definition
Network ACE (Agilent recommended)	ACE software is installed on a network node within the laboratory LAN infrastructure. Raw data locations are always captured in the equipment qualification report (EQR), which provides end to end traceability and a fully characterized data workflow in the delivery. This method requires collaboration with the customer to load ACE behind the customer firewall.
Local ACE	ACE software resides on an independent external drive that can be driven from the system controller, where the customer data system (CDS) resides. Because the external drive is connected to the CDS, the data integrity of this method is equivalent to that of the Network ACE delivery method. Raw data is imported directly into ACE by the Data Manager tool, with the data paths always captured in the report, which provides data traceability assurance.
Alternative (Requires pre-approval)	Pre-approval for this method is required to remove later questions on data integrity. ACE software is installed and run from a PC not directly connected to the CDS, such as the FSE laptop. System data files are transferred indirectly from the CDS to the laptop instead of directly like Network and Local ACE methods. NOTE: Software used in this method is qualified for data collection purposes; this method is <u>not</u> an option for software qualification.

EQR Storage

Select the checkbox below to authorize Agilent to store a copies of the EQRs generated by ACE for Agilent internal assessments. The intention of the assessment is to evaluate the delivery of the qualification service, with a focus to improve delivery and assess the appropriateness of data integrity measures. The storage is exclusively for the internal assessment by Agilent and is not shared with other organizations. It is not to be considered a backup for the EQR provided at qualification delivery.

Customer Approval of Alternative Method and EQR Storage

Authorize Agilent to use the alternative method (check for approval): ☐

Authorize Agilent to store EQRs for their internal assessment (check for approval): ☐

Approved By/Title:

Date Approved:

Comments:

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.
© Agilent Technologies, Inc. 2023
Published in USA



Introduction

With heightened scrutiny of data integrity, the Agilent Automated Compliance Engine (ACE) software must be able to access instrument-generated raw data files one of two ways: directly, using the connection between network nodes or with the server; and indirectly, through storage in a secure transfer location. (In this document, data integrity refers to the who, what, and where of data used in generating an ACE equipment qualification report, or EQR.)

ACE includes three main service delivery methods that address data integrity requirements; the rest of this document provides details to determine which one best fits a customer's needs.

Regardless of the delivery method, ACE features and delivery procedures are compatible.


Network ACE (Agilent Recommended)

Overview



ACE software is installed on a network node within the laboratory LAN infrastructure, which requires collaboration with the customer to load ACE behind their firewall. Raw data locations are always captured in the EQR, which provides end-to-end traceability and a fully characterized data workflow in the delivery.

Details



Typical Network ACE installation diagram

Installing ACE in a separate node (a.k.a. the host PC) on the same network as the system controller offers data traceability that is equivalent to an installation on the system controller itself. The system controller (where the CDS resides) and the ACE host PC are identified and seen by the server and subject to the customer's data access controls and general IT policies. The CDS's audit trail records data movements between nodes or between the client and server, and ACE's data traceability features identify the original data directory and therefore ensures end-to-end data traceability

The ACE host PC has a separate/partitioned drive for ACE software. During ACE's installation, two services are setup on the operating system (OS): one for security and the other as a watchdog. Because the ACE host PC sits on the network as a shared drive, engineers access ACE through the networked drive: ACE is not installed on ACE Virtual Viewer PCs.

Requirements

Installation

- Install on a host PC with a separate drive (different from that of the OS)
- Attach to a network that clients can access
- 500 GB
- NTFS format
- User has local administration rights
- Customer installation instruction document is available

Operational

- User has an ACE node logon with a minimum of power user rights permissions; user also has a personal ACE account and password added through the ACE licensing tool
- Up to 5 users with 3 open sessions each can access the NDA simultaneously
- Exception to ports 11121-11141 on ACE node, clients, and switch's/Smart Hubs to be open on the network

Local ACE

Overview



ACE software resides on an independent drive that can be driven from the system controller, where the CDS resides. Because the drive is connected to the CDS, this method's data integrity is equivalent to preferred 1 method's. Raw data is imported directly into ACE by ACE's Data Manager tool, and data paths are captured in reports to provide data traceability.

ACE software resides on an independent drive that can be driven from the system controller, where the CDS resides. Because the drive is connected to the CDS, this method's data integrity is equivalent to the Network ACE method. Raw data is imported directly into ACE by ACE's Data Manager tool, and data paths are captured in reports to provide data traceability.

Details

ACE is designed to run from a dedicated drive, without leaving a footprint on the host PC. Therefore, it can be connected directly to the system controller (where the CDS resides) without altering the system's qualification status. For additional protection, the drive can be driven by another host PC on the same network; also, the drive can remain on site with the customer for use by the Agilent Field Service Engineer (FSE) during service deliveries only.

Alternative Method

The ACE software is installed on and run from a PC not directly connected to the CDS, such as the FSE's laptop. System data files are transferred indirectly from the CDS to the laptop instead of directly like Network and Local ACE methods.

This method requires customer pre-approval to remove later questions on data integrity.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2023
Published in USA



Standard OQ Test Suite

This document describes the test program for qualifying ICP-OES instruments and their autosamplers. The ICP-OES tables list all tests run from software that controls the instrument; setpoints and limits cannot be changed and, for Agilent instruments, are internal to the controlling software.

ICP-OES: Agilent 5800, 5900, 5100, and 5110 Series

Note: Tests and test components are configuration specific as noted.

RV:	Radial Dual View
SVDV:	Synchronous Vertical Dual View
VDV:	Vertical Dual View

Test	Setpoint	Limits	
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data	
Preparation			
Plasma ignition	N/A	Within three attempts	
Detector calibration	N/A	Completes successfully	
Wavelength calibration	N/A	Completes successfully	
Instrument Test			
Subsystem Communication	N/A	Within acceptance criteria	
Air Flow	N/A	Within acceptance criteria	
Water Flow	N/A	Within acceptance criteria	
Gas Flows	N/A	Within acceptance criteria	
RF Generator	N/A	Within acceptance criteria	
Camera	N/A	Within acceptance criteria	
Optics	N/A	Within acceptance criteria	
		Models	
		5800 and 5900 Series	5110 and 5100 Series
Resolution	N (174.213 nm)	≤ 8.2	≤ 9.4
	As (188.980 nm)	≤ 8.2	≤ 8.2
	C (193.027 nm)	≤ 9.2	≤ 11.5
	Mo (202.032 nm)	≤ 8.2	≤ 8.2
	Cr (206.158 nm)	≤ 11.0	≤ 13.4
	Zn (213.857 nm)	≤ 8.7	≤ 8.7
	Pb (220.353 nm)	≤ 9.5	≤ 9.5
	Co (228.615 nm)	≤ 13.2	≤ 17.2
	Ba (230.424 nm)	≤ 9.4	≤ 9.4
	Mn (257.610 nm)	≤ 13.3	≤ 13.3
	Mn (260.568 nm)	≤ 15.0	≤ 20.3
	Cr (267.716 nm)	≤ 11.0	≤ 11.0
	Cu (324.754 nm)	≤ 24.0	≤ 25.0
	Cu (327.395 nm)	≤ 14.2	≤ 14.2
	Sr (338.071 nm)	≤ 25.0	≤ 33.5

Test	Setpoint	Limits			
Sensitivity, Axial (VDV, SVDV)	Ba (455.403 nm)	≤ 40.0		≤ 44.0	
	Sr (460.733 nm)	≤ 36.0		≤ 36.0	
	Ba (493.408 nm)	≤ 30.0		≤ 36.0	
	Ba (614.171 nm)	≤ 42.0		≤ 42.0	
	Ar (675.283 nm)	≤ 61.0		≤ 74.0	
	K (766.491 nm)	≤ 80.0		≤ 80.0	
	As (188.980 nm)	≥ 335.0		≥ 208.0	
	Se (196.026 nm)	≥ 205.0		≥ 159.0	
	Zn (206.200 nm)	≥ 2205.0		≥ 234.0	
	Zn (213.857 nm)	≥ 5325.0		≥ 1743.0	
	Cd (214.439 nm)	≥ 7465.0		≥ 4227.0	
	Pb (220.353 nm)	≥ 420.0		≥ 320.0	
	Mn (257.610 nm)	≥ 14625.0		≥ 10625.0	
	Cr (267.716 nm)	≥ 2210.0		≥ 1048.0	
	Cu (324.754 nm)	≥ 1620.0		≥ 19.0	
	Al (396.152 nm)	≥ 820.0		≥ 6.0	
	Ba (493.408 nm)	≥ 15080.0		≥ 60.0	
	K (766.491 nm)	≥ 5215.0		≥ 24.0	
Sensitivity, Radial	As (188.980 nm)	≥ 120.0 (VDV, SVDV)		≥ 46.0	
		≥ 110.0 (RV)			
	Se (196.026 nm)	≥ 60.0 (VDV, SVDV)		≥ 41.0	
		≥ 58.0 (RV)			
	Zn (206.200 nm)	≥ 660.0 (VDV, SVDV)		≥ 70.0	
		≥ 490.0 (RV)			
	Zn (213.857 nm)	≥ 2300.0		≥ 1421.0	
	Cd (214.439 nm)	≥ 2570.0 (VDV, SVDV)		≥ 522.0	
		≥ 1900.0 (RV)			
	Pb (220.353 nm)	≥ 130.0 (VDV, SVDV)		≥ 46.0	
		≥ 110.0 (RV)			
	Mn (257.610 nm)	≥ 5240.0 (VDV, SVDV)		≥ 3518.0	
		≥ 5600.0 (RV)			
	Cr (267.716 nm)	≥ 990.0 (VDV, SVDV)		≥ 379.0	
		≥ 1050.0 (RV)			
	Cu (324.754 nm)	≥ 560.0 (VDV, SVDV)		≥ 15.0	
		≥ 830.0 (RV)			
	Al (396.152 nm)	≥ 210.0 (VDV, SVDV)	≥ 3.4 (VDV, SVDV)		
		≥ 290.0 (RV)	≥ 4.6 (RV)		
	Ba (493.408 nm)	≥ 6180.0 (VDV, SVDV)	≥ 34.0 (VDV, SVDV)		
		≥ 9500.0 (RV)	≥ 55.0 (RV)		
	K (766.491 nm)	≥ 240.0 (VDV, SVDV)	≥ 1.8 (VDV, SVDV)		
		≥ 340.0 (RV)	≥ 3.0 (RV)		
Precision (RSD), Axial (VDV, SVDV)	As (188.980 nm)	≤ 1.5%			
	Se (196.026 nm)				
	Zn (206.200 nm)				
	Zn (213.857 nm)				
	Cd (214.439 nm)				
	Pb (220.353 nm)				
	Mn (257.610 nm)				
	Cr (267.716 nm)				
	Cu (324.754 nm)				

Test	Setpoint	Limits
Precision (RSD), Radial	Al (396.152 nm)	
	Ba (493.408 nm)	
	K (766.491 nm)	
	As (188.980 nm)	≤ 2.6%
	Se (196.026 nm)	≤ 1.5%
	Zn (206.200 nm)	
	Zn (213.857 nm)	
	Cd (214.439 nm)	
	Pb (220.353 nm)	≤ 2.6%
	Mn (257.610 nm)	≤ 1.5%
	Cr (267.716 nm)	
	Cu (324.754 nm)	
	Al (396.152 nm)	
	Ba (493.408 nm)	
	K (766.491 nm)	

ICP-OES: Agilent 700 Series

Note: Tests and test components are configuration specific as noted.

ACK:	Axial Cyclonic Spraychamber, K-style Nebulizer
ACO:	Axial Cyclonic Spraychamber, OneNeb Nebulizer
ACS:	Axial Cyclonic Spraychamber, Seaspray Nebulizer
ATO:	Axial Twister Spraychamber, OneNeb Nebulizer
ATS:	Axial Twister Spraychamber, Seaspray Nebulizer
RSMO:	Radial Sturman Masters Spraychamber, OneNeb Nebulizer
RSMV:	Radial Sturman Masters Spraychamber, V-groove Nebulizer
RTO:	Radial Twister Spraychamber, OneNeb Nebulizer
RTS:	Radial Twister Spraychamber, Seaspray Nebulizer

Test	Setpoint	Limits	
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data	
Preparation			
Plasma ignition	N/A	Within three attempts	
Detector calibration (71x-ES)	N/A	Completes successfully	
RF power check	N/A	Completes successfully	
Argon ratio check (72x-ES and 73x-ES)	N/A	Completes successfully	
Zinc wavelength position check	N/A	Completes successfully	
Dark current scan (72x-ES and 73x-ES)	N/A	Completes successfully	
Wavelength calibration	N/A	Completes successfully	
Hardware calibration (71x-ES)	N/A	Completes successfully	
Torch alignment	N/A	Completes successfully	
Instrument Test		Models	
		71x-ES	72x-ES
Resolution (72x-ES and 73x-ES)	Al (167.019 nm)	N/A	≤ 9.7
	N (174.213 nm)		≤ 9.4
	As (188.980 nm)	≤ 10.5	≤ 8.2
	C (193.027 nm)	N/A	≤ 11.5
	Mo (202.032 nm)		≤ 8.2
	Cr (206.158 nm)		≤ 13.4
	Zn (213.857 nm)	≤ 11.5	≤ 8.7

Test	Setpoint	Limits			
	Pb (220.353 nm)	N/A	≤ 9.5		
	Co (228.615 nm)		≤ 17.2		
	Ba (230.424 nm)		≤ 9.4		
	Mn (257.610 nm)		≤ 13.3		
	Mn (260.568 nm)		≤ 20.3		
	Cr (267.716 nm)		≤ 11.0		
	Cu (324.754 nm)		≤ 25.0		
	Cu (327.395 nm)		≤ 14.2		
	Sr (338.071 nm)	≤ 42.0	≤ 33.5		
	Ba (455.403 nm)	≤ 55.0	≤ 44.0		
	Sr (460.733 nm)	N/A	≤ 36.0		
	Ba (493.408 nm)		≤ 36.0		
	Ba (614.171 nm)		≤ 42.0		
	Ar (675.283 nm)	≤ 88.0	≤ 74.0		
	K (766.491 nm)	≤ 101.0	≤ 80.0		
		Models			
		710-ES		715-ES	
Signal Background Ratio	Pb (182.143 nm)	≥ 7.2		≥ 2.7	
	As (188.980 nm)	≥ 13.0		≥ 6.0	
	Se (196.026 nm)	≥ 7.5		≥ 3.8	
	Zn (206.200 nm)	≥ 60.0		≥ 11.0	
	Pb (220.353 nm)	≥ 12.0		≥ 3.4	
	Co (228.615 nm)	≥ 49.0		≥ 15.0	
	Ni (231.604 nm)	≥ 34.0		≥ 11.0	
	Cu (327.395 nm)	≥ 15.0		≥ 9.0	
	K (766.491 nm)	≥ 37.5		≥ 2.3	
		Models			
		720-ES 730-ES ACK, ACS, ACO	720-ES 730-ES ATO, ATS	725-ES 735-ES RSMO, RSMV	725-ES 735-ES RTO, RTS
Signal Background Ratio	Al (167.019 nm)	≥ 300	≥ 265	≥ 75.0	≥ 112.0
	Pb (182.143 nm)	≥ 11	≥ 9	≥ 6	≥ 8
	As (188.980 nm)	≥ 15.0	≥ 12.0	≥ 7.5	≥ 11.0
	As (193.696 nm)	≥ 7.5	≥ 6.0	≥ 6.0	≥ 9.0
	Se (196.026 nm)	≥ 9.0	≥ 7.0	≥ 4.5	≥ 6.5
	Zn (206.200 nm)	≥ 105	≥ 84	≥ 15	≥ 22
	Zn (213.857 nm)	≥ 190	≥ 150	≥ 110	≥ 150
	Cd (214.439 nm)	≥ 260	≥ 210	≥ 56.5	≥ 85.0
	Pb (220.353 nm)	≥ 15.0	≥ 12.0	≥ 3.8	≥ 6.0
	Mn (257.610 nm)	≥ 375	≥ 300	≥ 150	≥ 225
	Cr (267.716 nm)	≥ 52.5	≥ 45.0	≥ 22.5	≥ 34.0
	Cu (324.754 nm)	≥ 22.5	≥ 19.0	≥ 19.0	≥ 26.0
	Al (396.152 nm)	≥ 7.5	≥ 6.0	≥ 3.8	≥ 6.0
	Ba (493.408 nm)	≥ 75.0	≥ 60.0	≥ 75.0	≥ 112.0
	K (766.491 nm)	≥ 30.0	≥ 24.0	≥ 2.0	≥ 4.0
		All models			
QC Test - Accuracy	Same as precision setpoints	± 3.0%			
		Models			

Test	Setpoint	Limits				
		710-ES	715-ES	720-ES 730-ES	725-ES 735-ES	
QC Test - Precision	Al (167.019 nm)	N/A	N/A	≤ 1.50%	≤ 2.60%	
	Pb (182.143 nm)	≤ 1.60%	≤ 1.50%			
	As (188.980 nm)	≤ 1.00%	≤ 2.60%			
	As (193.696 nm)	N/A	N/A			
	Se (196.026 nm)	≤ 1.00%	≤ 2.60%			
	Zn (206.200 nm)	≤ 1.60%	≤ 1.50%			≤ 1.50%
	Zn (213.857 nm)	N/A	N/A			
	Cd (214.439 nm)					
	Pb (220.353 nm)	≤ 1.60%	≤ 1.50%		≤ 2.60%	
	Co (228.615 nm)			N/A	N/A	
	Ni (231.604 nm)					
	Mn (257.610 nm)	N/A	N/A	≤ 1.50%	≤ 1.50%	
	Cr (267.716 nm)					
	Cu (324.754 nm)					
	Cu (327.395 nm)	≤ 2.00%	≤ 1.50%	N/A	N/A	
	Al (396.152 nm)	N/A	N/A			
	Ba (493.408 nm)					
	K (766.491 nm)	≤ 1.00%	≤ 2.60%			
		Models				
		710-ES	715-ES	720-ES 730-ES	725-ES 735-ES RSMO RMSV	725-ES 735-ES RTO RTS
Detection Limits	Al (167.019 nm)	N/A	N/A	≤ 2.00	≤ 13.0	≤ 9.0
	As (188.980 nm)	≤ 10.0	≤ 65.0	≤ 10.0	≤ 65.0	≤ 45.0
	Se (196.026 nm)	N/A	N/A	≤ 13.0	≤ 80.0	≤ 50.0
	Mo (202.032 nm)			≤ 2.00	≤ 11.0	≤ 8.0
	Cd (214.439 nm)			≤ 0.50	≤ 5.50	≤ 4.0
	Pb (220.353 nm)	≤ 6.50	≤ 65.0	≤ 6.50	≤ 65.0	≤ 45.0
	Mn (257.610 nm)	≤ 0.20	≤ 0.90	≤ 0.20	≤ 0.90	≤ 0.60
	Cr (267.716 nm)	N/A	N/A	≤ 2.00	≤ 5.50	≤ 5.50
	Cu (324.754 nm)			≤ 2.00	≤ 5.50	≤ 5.50
	Al (396.152 nm)			≤ 6.50	≤ 13.0	≤ 13.0
	Ba (493.408 nm)			≤ 0.70	≤ 0.90	≤ 0.90
	K (766.491 nm)			≤ 20.0	≤ 265	≤ 265

ICP-OES: PerkinElmer Optima 4300DV, 5300DV, 7300DV, and 8000DV

Note: Test components are model specific as noted; **Scott** refers to an instrument with a Scott spray chamber and Gem tip crossflow nebulizer, as opposed to an instrument with a glass cyclonic spray chamber and concentric nebulizer.

Test	Setpoint	Limits
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data
Preparation		
OQ user creation	N/A	Completed
Software user paths defined	N/A	Completed
OQ support files loaded	N/A	Completed
Results data set and log files setup	N/A	Completed
Plasma ignition	N/A	Within three attempts
Detector calibration	N/A	Completes successfully

Test	Setpoint	Limits
Optics initialization (8000 DV only)	N/A	Completes successfully
Radial and axial torch alignment	N/A	Completes successfully
UV wavelength calibration (4300DV, 5300DV, 7300DV only)	N/A	Completes successfully
VIS wavelength calibration (4300DV, 5300DV, 7300DV only)	N/A	Completes successfully
Spectral Resolution		
8000DV	As (193.696 nm)	≤ 0.009
	Ni (231.604 nm)	≤ 0.011
	Ni (341.476 nm)	≤ 0.015
	Ba (455.403 nm)	≤ 0.020
4300DV, 5300DV, 7300DV DV	As (193.696 nm)	≤ 0.007
	Ni (231.604 nm)	≤ 0.008
	Ni (341.476 nm)	≤ 0.012
	La (408.672 nm)	≤ 0.020
	Ba (455.403 nm)	≤ 0.025
Precision (RSD)		
8000DV	Zn (206.200 nm)	≤ 1.0%
	Mg (280.271 nm)	
	Mg (285.213 nm)	
	Ba (455.403 nm)	
4300DV, 5300DV, 7300DV	As (193.696 nm)	≤ 1.0%
	Zn (213.856 nm)	
	Mn (257.610 nm)	
	La (379.478 nm)	
	Ba (455.403 nm)	
	Ba (493.408 nm)	
Axial Detection Limits		
8000DV	Tl (190.801 nm)	≤ 10 µg/L
	As (193.696 nm)	≤ 10 µg/L
	Se (196.026 nm)	≤ 5 µg/L
	Pb (220.353 nm)	≤ 3 µg/L
4300DV, 5300DV, 7300DV (*N/A for Scott)	Tl (190.800 nm)	≤ 10 µg/L
	As (193.696 nm)	≤ 10 µg/L
	Se (196.026 nm)*	≤ 5 µg/L
	Pb (220.353 nm)	≤ 3 µg/L
Radial Detection Limits		
8000DV	As (193.696 nm)	≤ 60 µg/L
	Zn (213.857 nm)	≤ 2 µg/L
	Mn (257.610 nm)	≤ 1 µg/L
	La (379.478 nm)	≤ 3 µg/L
	Ba (455.403 nm)	≤ 0.3 µg/L
	Ba (493.408 nm)	≤ 0.6 µg/L
4300DV, 5300DV, 7300DV	As (193.696 nm)	≤ 60 µg/L
	Zn (213.856 nm)	≤ 2 µg/L
	Mn (257.610 nm)	≤ 0.75 µg/L
	La (379.478 nm)	≤ 3 µg/L
	Ba (455.403 nm)	≤ 0.3 µg/L
	Ba (493.408 nm)	≤ 0.6 µg/L

Autosamplers and Switching Valve Accessories

This tests are qualitative; no setpoints or limits apply.

Test	Setpoint	Limits
Autosampler Operation		
Autosampler successfully moves to location(s) specified in software	N/A	Completed
Switching Valve Operation		
Automated test results are within acceptance criteria	N/A	Completed

Test Design and Rationale

Overview

Many GMP/GLP enforcement agency inspectors now ask firms to provide a risk assessment of their equipment and computer systems plus a science-based rationale for subsequent validation and qualification testing.

GENERAL RISK STATEMENT: Any laboratory chemical system used for raw material testing or final drug product / medical device testing in GMP or used in formal GLP studies will likely fall into a HIGH RISK category. This risk assessment will imply the need for IQ & OQ & on-going qualification. ANY USER SPECIFIC RISK ANALYSIS SUPERCEDES THIS GENERAL RISK STATEMENT.

The rest of this section outlines the science-based rationale for each test in the Agilent hardware OQ plus a brief test design and procedure description.

The recommended set of hardware OQ tests described in this EQP derives from Agilent's interpretation of FDA, USP, and GAMP guidelines and other authoritative expert literature.

OQ test design incorporates both modular and holistic testing, which is a proven and regulatory acceptable approach. When applicable, direct metrology is used to test pump flow rates and thermal-controlled column compartments, for example. Holistic chemical testing is used to evaluate critical instrument characteristics

When applicable, certified reference standards and calibrated equipment are used.

Considering the number of setpoints, parameters, and conditions of each recommended OQ test, the proven concepts of worst case, range, and representative have been applied. If a property or characteristic is known to have its worst performance at one end of a range of use, this is the setpoint that should be tested and other setpoints are not required. If a property or characteristic has no known worst case, testing at the high and low points of the range of use is required. If there are too many possible use cases and conditions to realistically test (and none is a worst case), a representative sample for test is the best approach.

All Systems

CDS Logon Verification

Description: For traceability, evidence of the logon used to collect qualification data must be provided.

Procedure: Service delivery engineer captures evidence that is automatically included with this test in the EQR.

ICP-OES: Agilent 5800, 5900, 5100, and 5110 Series

Preparation

Description: This preliminary test must be completed before the actual OQ tests. It verifies that:

- plasma ignited;
- detector was calibrated;
- instrument wavelength scale was calibrated.

Procedure: From the Agilent ICP Expert software (1) plasma is ignited and then extinguished; (2) while the plasma is off, the native software is used to calibrate the detector; (3) the instrument wavelength scale is calibrated.

Instrument Tests

Description: This test includes three automated instrument performance tests:

- Resolution
- Sensitivity
- Precision

Procedure: From the Agilent ICP Expert software, execute the Instrument Tests. The software prompts the user to aspirate the blank or standard solution as required. When the tests are completed, the software generates a report with the results.

ICP-OES: Agilent 700 Series

Preparation

Description: This preliminary test must be completed before the actual OQ tests. It verifies that:

- plasma ignited;
- (71x-ES only) detector was calibrated;
- RF power was checked;
- (72x-ES and 73x-ES only) argon ratio was checked;
- zinc wavelength position was checked;
- (72x-ES and 73x-ES) dark current was scanned;
- wavelength scale was calibrated;
- (71x-ES) hardware was calibrated;
- torch was aligned.

Procedure: From the Agilent ICP Expert II software (1) plasma is ignited; (2) all applicable functions are completed.

Instrument Tests

Description: This test includes three automated instrument performance tests:

- Resolution
- Signal Background Ratio
- QC Accuracy
- QC Precision
- Detection Limits

Procedure: From the Agilent ICP Expert II software, execute the Instrument Tests. The software prompts the user to aspirate the blank or standard solution as required. When the tests are completed, the software generates a report with the results.

ICP-OES: PerkinElmer Optima 4300DV, 5300DV, 7300DV, and 8000DV

Preparation

Description: This preliminary test must be completed before the actual OQ tests. It verifies that:

- prerequisites were completed (software's OQ user, user file paths, and OQ files have been installed; instrument and software log files are set up);
- plasma ignited;
- detector was calibrated;
- (8000DV only) optics initialization was completed;
- axial and radial torch views were aligned;
- (4300DV, 5300DV and 7300DV only) instrument UV wavelength scale was calibrated;
- (4300DV, 5300DV and 7300DV only) instrument VIS wavelength scale was calibrated.

Procedure: From the instrument's software (1) prerequisites are completed; (2) plasma is ignited; (3) detector is calibrated; (3) 8000DV instrument's optics are initialized; (4) instrument's axial and radial torch views are aligned; (5) 4300DV, 5300DV, and 7300DV instrument's UV and VIS wavelength scales are calibrated.

Spectral Resolution

Description: This test measures the resolution, or full width at half height (FWHH), of the relevant element's spectral peak.

Procedure: A standard solution containing the relevant elements is aspirated using a dedicated method. The software is set up to log the three replicate resolution measurements for each element. The highest resolution measurement is recorded in ACE as the result.

Precision

Description: This test determines the precision (% RSD) of the relevant element's emission signal.

Procedure: A standard solution containing relevant elements is aspirated using a dedicated method. The method measures a number of replicates and calculates the % RSD for each element. The calculated result is then recorded in ACE.

Axial Detection Limits

Description: This test determines the minimum concentration at which the relevant element can be detected in the axial viewing mode.

Procedure: A dedicated method is used to aspirate a 2% nitric acid blank and then a multi-element calibration standard solution to formulate a calibration curve. The sample introduction system is then flushed with a 10% nitric acid wash solution. The 2% nitric acid blank solution is aspirated again and the replicate measurements are interpolated from the calibration curve. Statistical calculations are then applied to the interpolated results and the calculated results are recorded in ACE.

Radial Detection Limits

Description: This test determines the minimum concentration at which the relevant element can be detected in the radial viewing mode.

Procedure: A dedicated method is used to aspirate a 2% nitric acid blank and then a multi-element calibration standard solution to formulate a calibration curve. The sample introduction system is then flushed with a 10% nitric acid wash solution. The 2% nitric acid blank solution is aspirated again and the replicate measurements are interpolated from the calibration curve. Statistical calculations are then applied to the interpolated results and the calculated results are recorded in ACE.

Autosamplers and Switching Valve Accessories

Autosampler Operation

Description: This test verifies the functional operation of the installed autosampler.

Procedure: From the instrument software, the autosampler probe is directed to several sample locations and a pass or fail is recorded in ACE.

Switching Valve Operation

Description: This test verifies the functional operation of the installed switching valve.

Procedure: From the instrument software, run the automated switching valve test and follow the prompts. When the tests are completed, the software generates a report with the results and a pass or fail is recorded in ACE.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2023
Published in USA

FYdcfhUbX'8Y]j YfmCdHcbg

V@ÁaóÁ[[, Ág &`á^•Á] ǻ }•Á^|&óáÁ!Ág &`•ǻ }Á@ÁÖÜÛËÓ Á^æ|ǻ^] ^æ^áÁ~} Á[*•Ág áÁg}•æǻ } Á[*•Ág^Ág, æ•
ǻ &`á^áÁg|Á c@!Á] ǻ }•Ág Á^Ág &`á^áÁ!Ág &`á^áÉ

ǻ|!Á@] @} æ^áÁ~•c{ Ác] ^•Á} |`DÁáá^!^} ǻ] ǻ }•Ág^Ág•^} Á|!Á@!ǻ æ^ Ág áÁ~]] [|!ááÁ~•c{ Ác] ^•Ág^Á!ǻ æ^ Á^•c{
[] ǻ }•Ág^Á^•áÁ! Á|!Á c@Á&@ǻ~^•É

ÄÜ@, Ág@[{ æ[*!æ•
ÄÜ@, Á@æ^!Ág áÁ[c!Á} Ág ç!
ÄÜg &`á^Á^] ^æ^áÁ~} Á[*•
ÄÜg &`á^Á!g}•æǻ } Á[*•

GY'WYX'G][bUhi fY'CdHcbg

Üæ•KÖÜÚÁ Á| ǻ[&^á
ÄÜ^ [|!ǻ * Ágæg &^Ág Á[[, ^áÁg Ág ÁÖÜÜ

7 i ghca Yf`5 ddfcj U

Pađ ^K
Văđ^K
Öæ^K
Ùă}æ`!^K

Pađ ^K
Văđ^K
Öæ^K
Ùă}æ`!^K

Pađ ^K
Văđ^K
Öæ^K
Ùă}æ`!^K

Pađ ^K
Văđ^K
Öæ^K
Ùă}æ`!^K

